
Excavators


SPOONS: For Removal of Carious Dentin



MEC-101-E1

Pluggers / Condensers


Used to Pack Filling Materials into Cavity Preparations




MPL-110- 1/2,ANDREW

CARVERS

Used to Crave anatomical features and trim excess materials.




Hollow Handls (A)




Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Burnishers


Designed to condense, smooth, carve and polish amalgam.
Gold Titanium Coated Points.



MBU-101-18

PK-Thomas

Used for waxing procedures and techniques.



MPKT-101 **MPKT-102** **MPKT-103** **MPKT-103R** **MPKT-104** **MPKT-105**

Carving Knives

Used for trimming excess filling material, flash and overhangs.



MCK-113 / 20 Esthetic
For anterior teeth. Sharp, thin blade allows access to all surfaces. The offset angle provides universal adaptability



MCK-114 / 21 Esthetic
For posterior teeth. Thin, sharp offset angle provides access to surfaces where linear finishing strips would not be effective




MCK-115 / 6 Tanner
For posterior teeth. Sharp offset angle provides access to many surfaces



MCK-116 / 12 Novatech Cement Remover
Combines a sickle-shaped scaler with a flat blade for removal of excess resin. provides access to many surfaces cement. or porcelain flasn
The narrow chisel removes excess interproximal material with a push stroke


Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Composite Filling Instruments



MCG-101 /#1 Small Universal



MCG-102
#2 Medium Universal



MCG-103
#4 Small Universal



MCG-104
#3 Large Universal



MCG-105
#5 Medium Universal



MCG-106-Boghosian
Combination of Thin, Knife-Shaped Blade



MCG-107-8A
Use for packing gingival retraction cord.



MCG-108-Boghosian
Used for measuring composite layers.



MCG-109-Contouring
Used for shaping of inclines, planes
for Anterior and posterior restorations.



MCG-110
Quétin



MCG-111-8A
#2 Round Arrow

Composite Filling Instruments Set



MCB-125



MCB-126



MCBS1-01
Set of 6 in Cassette



MCB-127



MCB-128




MCB-129



MCB-130

Amalgam Carrier



MS15-110
Mini/Regular 90°


Amalgam Well



MSB-10D



MS15-111
Jumbo/Regular 30°
20 cm



MS15-112
Large/Regular 45°
20cm



MS15-113
mm,2,00 ø
15cm



MS15-114
mm,2,00 ø
15cm



MS15-116
mm,2,00 ø
15cm



MS15-116