Hard Wire Cutters



MS36-101
Hard wire Cutter



MS36-102
Hard wire Cutter



MS36-103
Hard wire Cutter



MS36-104
Hard wire Cutter


WIRE CUTTER TC GOLD

Tungsten Carbide (Tip Insert)



MS36-01-16 cm
NS36-02-18 cm
wire Cutter
TC Gold
Double Action

Hand Instruments



Ortho light scaler double ended — **MS36-105**
Twiel-on elastic module instrument — **MS36-106**
Ortho Light double ended ligature director — **MS36-107**
Ortho Light double ended ligature scaler/director tivate coil spring are remove modules with one instruments — **MS36-108**
Schure Band Seater — **MS36-109**
Band seater auto cleavable — **MS36-110**
Mershon Band Seater — **MS36-111**

Uses of Articulator

Primary purpose:

1) To hold opposing casts in a predetermined fixed relationship
2) To open and close. 3) To produce border and inter- border diagnostic sliding motions of the teeth similar to these in the mouth.







MSAR-1001
Technical Datas:
Intercondylar distance 120 mm. / condylar inclination fixed at 30°
Bennet angle fixed at 7° / Incisal table inclination 10°
Distance between mounting plates 85 mm. / W eight 400 gr.
Opening 180°

MSAR-1002



Black

MSAR-1003



Blue

MSAR-1004



Yellow

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co