
Excavators


SPOONS: For Removal of Carious Dentin



MEC-101-E1

Pluggers / Condensers


Used to Pack Filling Materials into Cavity Preparations




MPL-110- 1/2,ANDREW

CARVERS

Used to Crave anatomical features and trim excess materials.




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)



Silicon Handls (E)





MCP-116 / BN1
Thin blade and rounded contour facilitates use in both thick and thin tissues without catching or dropping cord. Bilateral notch allows placement in limited access areas.



MCP-117 / 1 Yardley
Round non-serrated working end.



Serated **MCP-118** Non-Serated
CSI-1



MCP-119
7 Guyer



Serated **MCP-120** Non-Serated
113

Gingival Retractors

Protects tissue during cavity preparations and placement of composite restorations




MGR-120 / 1 Meinershagen / For mandibular inscisors



MGR-121 / 2 Meinershagen
For maxillary and mandibular premolars and canines



MGR-122 / 3 Meinershagen
For maxillary central incisors and wide canines



MGR-123 / 4 Meinershagen
For all molars.



MGR-124
10 Goldman-Fox


Composite Filling Instruments

Also Available in Titanium Coated Colored Points.

Vertical and horizontal combination.




MCF-101



MCF-102

MCF-103




MCF-104

MCF-105



MCF-106

MCF-107



MCF-108

MCF-109




MCF-110
Woodson




MCF-111
W1

Vertical in and offset angled.



MCF-112


MCF-113



MCF-114

MCF-115


Hollow Handls (A)



Titanium Coated Handls (B)



Solid Handls (D)



Silicon Handls (E)

