
Excavators


SPOONS: For Removal of Carious Dentin



MEC-101-E1

Pluggers / Condensers


Used to Pack Filling Materials into Cavity Preparations




MPL-110- 1/2,ANDREW

CARVERS

Used to Crave anatomical features and trim excess materials.




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)



Silicon Handls (E)





MCP-116 / BN1
Thin blade and rounded contour facilitates use in both thick and thin tissues without catching or dropping cord. Bilateral notch allows placement in limited access areas.



MCP-117 / 1 Yardley
Round non-serrated working end.



Serated **MCP-118** Non-Serated
CSI-1



MCP-119
7 Guyer



Serated **MCP-120** Non-Serated
113

Gingival Retractors

Protects tissue during cavity preparations and placement of composite restorations




MGR-120 / 1 Meinershagen / For mandibular inscisors



MGR-121 / 2 Meinershagen
For maxillary and mandibular premolars and canines



MGR-122 / 3 Meinershagen
For maxillary central incisors and wide canines



MGR-123 / 4 Meinershagen
For all molars.



MGR-124
10 Goldman-Fox


Composite Filling Instruments

Also Available in Titanium Coated Colored Points.

Vertical and horizontal combination.




MCF-101



MCF-102

MCF-103




MCF-104

MCF-105



MCF-106

MCF-107



MCF-108

MCF-109




MCF-110
Woodson




MCF-111
W1

Vertical in and offset angled.



MCF-112


MCF-113



MCF-114

MCF-115


Hollow Handls (A)



Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Composite Filling Instruments



MCG-101 /#1 Small Universal



MCG-102
#2 Medium Universal



MCG-103
#4 Small Universal



MCG-104
#3 Large Universal



MCG-105
#5 Medium Universal



MCG-106-Boghosian
Combination of Thin, Knife-Shaped Blade



MCG-107-8A
Use for packing gingival retraction cord.



MCG-108-Boghosian
Used for measuring composite layers.



MCG-109-Contouring
Used for shaping of inclines, planes
for Anterior and posterior restorations.



MCG-110
Quétin



MCG-111-8A
#2 Round Arrow

Composite Filling Instruments Set



MCB-125



MCB-126



MCBS1-01
Set of 6 in Cassette



MCB-127



MCB-128




MCB-129



MCB-130

Amalgam Carrier



MS15-110
Mini/Regular 90°


Amalgam Well



MSB-10D



MS15-111
Jumbo/Regular 30°
20 cm



MS15-112
Large/Regular 45°
20cm



MS15-113
mm,2,00 ø
15cm



MS15-114
mm,2,00 ø
15cm



MS15-116
mm,2,00 ø
15cm



MS15-116