MRSHAHco

QUALITY DENTAL INSTRUMENTS

A Better Life Start With a Beautiful Smile







For Inquiries

+49 174 8102485

info@mrshah.co

www.mrshah.co

Mrshah Co. (GmbH)

DENTAL INSTRUMENTS

Catalouge

A Little About Us:

Mrshah Co. Offers A Wide Variety Of General Dental Instruments. Our High Quality Instruments Provide Dentists With The Proper Instrumentation For A Variety Of Dental Procedures.

We Please To Present The New Edition Dental Instruments Catalogue With Hope That It Would Be An Excellent Informative Source To Our Esteemed Customer, Distributors And Dental Professionals Worldwide Which Includes Following Sections,Diagnostics, Periodontology, Surgery, Implant, Extraction, Restorative & Orthodontics.

Mrshah Co. Is ISO-9001-2008, ISO-13485-2012,CE, & cGMP (FDA-USA) Certified Company And Head Office Located In Frankfurt Am Main, Germany And Manufacturing Unit In Small Industrial Estate Of Sialkot,Pakistan . Mrshah Co. Is A Leading Manufacturer And Reliable Source Of Excellent Quality Dental Instruments.

Due To Continuous Policy Of Innovation And Further Development In The Field. Our Products Are Better Designed, More Comfortable And Long-Lasting.
We At Mrshah Co. Successfully Running

Family Owned Business For Last 26years With The Support Of Great Team, Whose Skilled Techniques AndTraditional Approach Towards Modern Surgery,Make All Possibilities To Acclaim Us Quality Instruments Manufacturer.

We Are Also Very Grateful For The Loyal Support Given By Our Customers/Distributors Who Have Been Selling Our Instruments In Their Market With Pride And Honor.

We Are Pledge To Serve More In The Years Ahead To Our Quality Minded Customers With Much More Higher Standard Of Instruments And Services.

We Offer A Comprehensive Portfolio Of Reusable And Disposable Instruments For Dental Surgeries, Dental Laboratories, And For Dental Students.

We Also Welcome Your Enquires For Instruments Not Listed In This Catalogue Both In Surgical & Dental Instruments.

At Mrshah Co. We Feel That Good Communication Is Key To Any Successfully Business Relationship. PleaseFeel Free To Contact Us By The Methods That Are TheMost Convenient For Your Business And We Will Do OurBest To Always Provide Great Quality Products AndServices.

Yours Sincerely.
Mrshah Co (GmbH) Management Team.

MRSHAH CO (GmbH)
WARRANTY

Mrshah Co warranties all its products against poor craftsmanship and materials. Any product found defective will be repaired or replaced free of charge. This warranty does not cover any products which become defective through misuse, negligence and mishandling.
MRSHAH CO (GmbH)



Index

01 Extraction Forceps (English | American | Children)  Page | 01-09

02 Bone Rongeurs  Page | 10

03 Root Elevators  Page | 11-16

04 Dental Anesthesia Syringes  Page | 17

05 Mirrors Handles & Retainer Matrix  Page | 18

06 Dental Tweezers  Page | 19

07 Dental Scissors  Page | 20

08 Needle Holders  Page | 21-22

09 Explorers & Probes  Page | 23

10 Excavators / Pluggers / Condensers / Carvers  Page | 24

11 Burnishers/ Hollenback / Carving Knives  Page | 25

12 Restorative  Page | 26-27

13 Dental Scalers & Curettes  Page | 28-30

14 Periodontal Knives  Page | 31

15 Dental General Instruments  Page | 32


16 Rubber Dam Instruments  Page | 33

17 Dental Crown Instruments  Page | 34

18 Dental Implants Instruments  Page | 35-37

19 Orthodontics Instruments  Page | 38-40


20 Impression Trays  Page | 41

21 Dental Hollow Wares / Cassettes  Page | 42-43






CHILDREN& ADULT TEETH



Permanent Teeth




Requirment of the Dental Forceps

INTRA-ALVEOLAR EXTRACTION

Dental extraction forceps has three components which is blades and handles united together by a hinge joints.
1-The serrated handle must possess a suitable size to rest comfortably in the operator's hand and should long enough to afford use to strong and steady extraction movements, Also it must be suitable in shape and designs to suit the area of the tooth.







Extracting Forceps (Set of 6 Pcs)

MRSHAHco
QUALITY DENTAL INSTRUMENTS

Set Of 6

Cross serrations for a non-slip grasp. Easily enter the sub gingival.


Unique Desgin.

Universal Atraumatic Extraction Forceps Set Premium (Adult / Pediatric / Children)



MSS1-01
Set Of 6 Pcs





MSA1-154 **MSA1-155** **MSA1-156** **MSA1-157** **MSA1-158**

MSA1-153 Premium Nature Atraumatic Universal Extraction Forceps Lower Anterior
MSA1-154 Premium Nature Atraumatic Universal Extraction Forceps Upper Anterior
MSA1-155 Premium Nature Atraumatic Universal Extraction Forceps Universal Upper Premolars
MSA1-156 Premium Nature Atraumatic Universal Extraction Forceps Universal Lower Premolars
MSA1-157 Premium Nature Atraumatic Universal Extraction Forceps Lower Molars
MSA1-158 Premium Nature Atraumatic Universal Extraction Forceps Upper Molars



MS1-101
Fig.1
Upper Centrals
and Canines



MS1-102
Fig.2
Upper Leternals
and Canines



MS1-103
Fig.3
Upper Incisors
and Canines



MS1-104
Fig.4
Lower Incisors
and Canines



MS1-105
Fig.7
Upper Premolars



MS1-106
Fig.8
Lower Premolars



MS1-107
Fig.13
Lower Premolars



MS1-108
Fig.17
Upper Molars
Right



MS1-109
Fig.18
Upper Molars
Left



MS1-110
Fig.19
Upper Third
Molars



MS1-111
Fig.20
Lower Third
Molars



MS1-112
Fig.21
Lower Molars



MS1-113
Fig.22
Lower Molars



MS1-114
Fig.30
Upper Roots



MS1-115
#33-Lower Roots



MS1-116
#33A-Lower Roots



MS1-117
#73A-Lower Molars



MS1-118
#79-Lower Third Molars



MS1-119
#45-Lower Roots



MS1-120
#44-Upper Roots



MS1-121
#51A-Upper Roots



MS1-122
#67A-Upper Third Molars



MS1-125
#73-Lower Molars



MS1-126
#74-Lower Roots



MS1-127
#74N-Lower Roots



MS1-123
Routurier (#22 ½ L)
Lower Molars and
wisdoms, left

WISDOM EXTRACTING FORCEPS




MS1-124
Routurier (#22 ½ R)
Lower Molars and
wisdoms, right



MS1-128
#86A-Lower Molars

DEEP GRIPPING FORCEPS

Parallel Gripping Beak
Aggressive edge serrations
Long axis serrations within the beaks
Thick beaks for superior penetration
Contact the root remaning coronal structure
Used on partially intact, broken down or sectioned teeth



MRSHAHco
QUALITY DENTAL INSTRUMENTS

ORAL SURGERY





No Fracture | No Tissue Damage | Extraction with Easy | Slim Light Weight Pattern



MS1-129
#35
Upper Permolars
Deep-gripping

Thin Tip



MS1-130
#34
Upper Incisors and
Premolars
Deep-gripping

Wide Tip



MS1-131
#34A
Upper Centrals / Cains
Deep-gripping

Wide Tip



MS1-132
#35A
Upper Incisors and
Premolars
Deep-gripping

Thin Tip



MS1-133
#36
Lower Premolars
Deep-gripping

Wide Tip



MS1-134
#36A
Lower Permolars
Deep-gripping




MS1-135
#336
Lower Permolars
Deep-gripping



MS1-136
F-7C
Lower Molar




MS1-137
#9
Apex Foeceps
Lower Universal
Deep Root Pick



MS1-138
#79N
Lower Third Molars
Deep-Gripping



MS1-139
#79N
Lower Third Molars



MS1-140
#67
Upper Third Molars



MS1-141
#67A
Upper Third Molars



MS1-142
#67N
Upper Third Molars
Deep-gripping

Extracting Forceps (Set of 10 Pcs)



MS1-143 Fig.1, Upper Leterals and Canines



MS1-144
Fig.7
Upper Premolars



MS1-145
Fig.13
Lower Premolars



MS1-146
Fig.17 Upper Molars Right



MS1-147
Fig.18 Upper Molars Left



MS1-148
Fig.22
Lower Roots



MS1-149
Fig.51A
Lower Roots



MS1-150
Fig.67 Lower Third Molars



MS1-151
Fig.74N
Upper Molars



MS1-152
Fig.79 Upper Third Molars



MSS1-02- Set of 10 Pieces
Fig.1, Fig.7, Fig.13, Fig.17
Fig.18, Fig.22, Fig.51A, Fig.67
Fig.74N, Fig.79

Universal Pediatric Forceps



MS1-153
Fig.365/150SK
Upper anteriors and Roots



MS1-154
Fig.360/17SK
Lower Molars



MS1-155
Fig.370/151SK
Lower anteriors and Roots

Root Pick Extracting Forceps



MS1-156
Forceps Fig.301



MS1-157
Forceps Fig.300

Perforated Light Weight Handles with Strong Serrated Gripping Action.



MS2-101
Fig.1
Upper Anterior

Side View



MSS2-01
Set of 10 Pieces
Fig.1, Fig.23, Fig.174
Fig.88L, Fig.88R, Fig.74
Fig.106, Fig. 174, Fig.859

MS2-102
Lower Molar
Fig.23 Cow Horn



MS2-103
Fig.35
Upper Premolar



MS2-104
Fig.88L
Upper Molar Left



MS2-105
Fig.88R
Upper Molar Left



MS2-106
Fig.74
Lower Molar



MS2-107
Fig.106
Lower Incisors



MS2-108
Fig.174
Extra Narrow
Lower Anterior
Apical Retention
Forceps



MS2-109
Fig.74
Upper Molars



MS2-110
Lower Roots



BABY EXTRACTING FORCEPS


To avoid fear of small patients, **MR SHAH** forceps are that small tofit completely inside a hand. For an easy use they are additionally equipped with a spring.**MR SHAH** extracting forceps assure a firm, accurate grasp of primary tooth and root surface for smooth, seamless extraction reducing the risk of complicattions. Hand crafted by Surgical Grade Non-Rusty Stainless steel superior perfomance and longevity.

Think Teeth

Keep up good dental habits at every age.



Pediatric Tooth Extracting Forceps Sets



MS3-101
English Pattren Klein #137 Pedodontic, Upper Incisors Spring handle



MS3-102
English Pattren Klein #5 Pedodontic Lower Incisors.



MS3-103
English Pattren Klein #7 Pedodontic Lower Roots,



MS3-104
English Pattren Klein #6 Pedodontic Lower Molars, Spring Handle



MS3-105
English Pattren Klein #139 Pedodontic Upper Premolars



MS3-106
English Pattren Klein #51S Pedodontic, Upper root




MS3-107
English Pattren Klein #3 Pedodontic, Upper Molars Universal



MSS3-01
Set of Pieces

Pediatric Tooth Extracting Forceps Sets



MS3-108
Englisg Pattren Klein #137 Pedodontic, Upper Incisors Spring handle



MS3-109
English Pattren Klein #5 Pedodontic, Lower Incisors.



MS3-110
English Pattren Klein #7 Pedodontic, Lower Roots,



MS3-111
English Pattren Klein #6 Pedodontic, Lower Molars, Spring Handle



MS3-112
English Pattren Klein #139 Pedodontic, Upper Premolars



MS3-113
Englisg Pattren Klein #51S Pedodontic, Upper root



MS3-114
Englisg Pattren Klein #3 Pedodontic, Upper Molars Universal




MSS3-02
Set of 7 Pieces





MS4-101
Fig.62
incisor, premolars,
decidous teeth



MS4-102
#10S
upper roots



MS4-103
Harris
Fig.18R



MS4-104
Harris
Fig.18L



MS4-105
Fig.53R
upper molars, right



MS4-106
Fig.53L
upper molars, left



MS4-107
NEVIUS
Fig.88R
upper molars, right



MS4-108
Fig.24
upper molars



MS4-109
#210S
upper third molars



MS4-110
#99C Kells
upper incisors,
canines,premolars



MS4-111
Nevius
Fig.88L
upper molars, left



MS4-112
Tomes Fig.69
upper and lower
roots



MS4-113
Hull Fig. 101
Deciduous teeth



MS4-114
Fig. 150 AS



MS4-115
Fig. 222
Lower wisdoms



MS4-116
Fig. 17
Lower molars



MS4-117
Fig. 16
Harris Lower molars,
COW HORN beak



MS4-118
Fig. 222
Lower wisdoms



MS4-119
Fig. 69
Tomes Fragments,
upper and lowers



MS4-120
Fig. 286
Upper roots,
premolars and incisors



MS4-121
Fig. 171
Wisdom and molar roots



MS4-122
Fig. 166
Universal Application



MS4-123
Fig. 169
Centrals laterals premolars
lower canines and all roots



MS4-124
Fig. 24
upper molars



MS4-125
Fig. 151
Cryer Lower Premolars,
incisors, roots, universal



MS4-126
Fig. 76
Upper roots



MS4-127
Tweezer Albrecht



MS4-128
Scissor Handle
Lower Wisdoms

MRSHAHco
QUALITY DENTAL INSTRUMENTS



BONE RONGEUR

INSTRUMENTS FOR CUTTING BONES

- Most commonly used
- Sharp blades- squeezed together cutting/pinching through bone
- Leaf spring between the handle : instrument opens when hand pressure is released
- Repeated cuts without manually reopening



MS5-101
Blumenthal 90°
15cm



MS5-102
Blumenthal 45°
15cm



MS5-103
Blumenthal 30°
15cm



MS5-104
Friedman
14cm



MS5-105
Mini Friedman
14.5cm



MS5-106
Mini Friedman 90°
14.5cm

SOFT TISSUE NIPPERS



MS5-113 - 10cm
Compact, sharp nippers with smooth action useful for accentuating interproximal gingival contour during soft tissue surgery.

DOUBLE ACTION RONGEURS



MS5-107
Double Action
Side-Cutting
17 cm



Boehler
Double Action
STR
MS5-108 - 15 cm-2mm
MS5-109 - 15 cm-3mm
CVD
MS5-110 - 15 cm-2mm
MS5-111 - 15 cm-3mm



MS5-112
Beyer
Double Action
18 cm

KERISON RONGEUR




Quickly and Easily Enlarge
The Sinus Window For
Increased Access
& Visibility

MSC5-150-2mm
MSC5-151-3mm
Kerison Rongeur

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

HANDLES FOR ROOT ELEVATORS



Elevators Classification

- I-According to use:

1-Elevators designed to remove the entire tooth, straight elevators, hospital pattern and winter elevator

2-Elevators designed to remove roots broken off at the gingival line e.g. Apexo elevator, Coupland and lido lavien elevators

3-Elevatos designed to remove roots broken off half way to the apex e.g. curved elevator hospital patters, winter elevator and Apexo elevator



HANDLES FOR ROOT ELEVATORS POINTS / A-B-C-D-E-F-G-H


A
Hollow Handle (Dull Finish)


B
Hollow Handle (Short Blast Coated)


C
Hollow Handle


D
Hollow Handle (Sand Blast)


E
Hollow Handle


F
Hollow Handle


G
Hollow Handle


H
Hollow Handle

Flexible Periotomes

Gold Titanium Coated Points with Micro Serration



MS6-101
Flexible Periotome 2mm CVD

MS6-102
Flexible Periotome 2mm STR

Root Screw Extractor



MS6-103

MR SHAH Flexitome (Set of 6 Pcs)

Super Flexibility, Micro Serrated, Special Alloy, Hard To Be Broken



MS6-104
Straight 3mm



MS6-105
Straight 2mm

MS6-106
Curved 2mm

MS6-107
Curved 2mm

MS6-108
Side Cutting 2mm

MS6-109
Back Cutting 2mm

Flexitome Cassette



MSS6-01
Flex-PDL Set of 6

Apexo Root Elevators



SM6-110-Apexo
2.5mm-Straight



SM6-111
Apexo
2.7mm
Straight



SM6-112
Apexo
3.0mm
Straight



SM6-113
Apexo
3.5mm
Straight



SM6-114
Apexo
2.5mm
Curved



SM6-115
Apexo
2.7mm
Curved



SM6-116
Apexo
3.0mm
Curved

Bein Root Elevators



SM6-117-Bein
2.0mm

2.5mm



SM6-118
Bein 2.5mm

3.0mm



SM6-119
Bein 3.0mm

4.0mm



SM6-120
Bein 4.0mm

Seldin Root Elevators

3.5mm



SM6-121
Seldin1
3.5mm

4.0mm



SM6-122
Seldin34S
4.0mm

5.0mm



SM6-123
Seldin34
5.0mm

Apical Root Elevators



SM6-124
Fig.1
Straight



SM6-125
Fig.4
Left



SM6-126
Fig.5
Right

Cryer Root Elevators



SM6-127
Cryer Left
Small



SM6-128
Cryer Left
Small



SM6-129
Cryer L
Medium



SM6-130
Cryer R
Medium



SM6-131
Cryer L
Large



SM6-132
Cryer R
Large

Coupland Root Elevators



SM6-133
Coupland
2.5mm

3mm



SM6-134
Coupland
No.1

3.5mm



SM6-135
Coupland
No.2

4mm



SM6-136
Coupland
No.3

5mm



SM6-137
Coupland
No.4

Crane Pick



SM6-138
E8

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

Lindo Root Elevators

MICRO SERRATED



SM6-139



SM6-140
25
LINDO LEVIAN
2.5mm
Straight



SM6-141
35
LINDO LEVIAN
3.0mm
Straight



SM6-142
45
LINDO LEVIAN
4.0mm
Straight



SM6-143
1
LINDO LEVIAN
2.5mm
Curved



SM6-144
2
LINDO LEVIAN
3.5mm
Curved



SM6-145
3
LINDO LEVIAN
4.0mm
Curved

Bernard Root Elevators



SM6-146
Bernard Straight



SM6-147
Bernard
Curved



SM6-148
Bernard Serrated
in Offset



SM6-149
Bernard Serrated
out Offset



SM6-150
Bernard
Left



SM6-151
Bernard
Right

Claw Lift Root Elevators




SM6-154
Claw-Lift
Left

SM6-155
Claw-Lift
Right



SM6-153
Winter
11R



SM6-152
Winter
11L

Heidbrink Root Tip Picks



SM6-156
Heidbrink
Left



SM6-157
Heidbrink
Right



SM6-158
DEL H1



SM6-159
DEL H2



SM6-160
DEL H3

MICRO SERRATED



SM6-161
Heidbrink
Straight



SM6-162
Heidbrink
Left



SM6-163
Heidbrink
Right




SM6-164
Heidbrink
Double
Ended

For those clinic-ions who are interested the entire luxating Elevators line. **MR SHAH** has created set with seven micro serrated instruments for Automatic extraction

PDL Elevators (Set of 7 Pcs)



PDL Elevators (Set of 7 Pcs)



PDL Elevators (Set of 7 Pcs)




Bernard Flexible (Set of 4 Pcs)



TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

Luxation Elevators



MS6-190 **MS6-191** **MS6-192** **MS6-193** **MS6-194**



MS6-195 **MS6-196**




MSS6-05 Luxation set of 7 Pcs

Luxation Elevators



MS6-197
2.0 mm
Straight



MS6-198
3.0 mm
Straight

MS6-199
4.0 mm
Straight

MS6-200
5.0 mm
Straight

MS6-201
2.0 mm
Curved

MS6--202
3.0 mm
Curved

MS6-203
4.0 mm
Curved

MS6-204
5.0 mm
Curved

Potts Luxation Elevators



MS6-205
Left

MS6-206
Right

Surgical Curettes



MS6-207
3 with
hollow handle

Syndesmotorme Chompret



MRC-101 **MRC-102** **MRC-103** **MRC-104**

Lucas Surgical Curettes

MICRO SERRATED



MRC-105
Lucas
3.0 mm

MRC-106
Lucas
3.0 mm

MRC-107
Lucas
3.0 mm

Luxation Set, Delicated Points



MS6-208
straight blade
2mm



MS6-209
curved blade
3mm



MS6-210
straight blade
3mm



MS6-211
Contra angle
3mm



MS6-212
Inverted curved
3mm



MS6-213
Straight blade
5mm



MS6-214
curved blade
5mm



MS6-215
Straight dual edge
3mm/3mm



MS6-216
Straight dual edge
3mm/1.5mm



MS6-217
Right to the
bending 3mm



MS6-218
Left to bend
3mm



MS6-219
Right to the
bending 5mm



MS6-220
Left to bend
5mm

Ergonomic Handles Elevators Set



MS6-221
Apical Elevator
left



MS6-222
Apical Elevator
Right



MS6-223
Straight 3mm



MS6-224
Cryer Elevator
Left



MS6-225
Cryer Elevator
Right



MS6-226



MS6-227



MS6-228



MS6-229




MSS6-06

Root Fragment Removal elevators, Set of 4 (Twisted Style)

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

Self Aspirating Local Anesthesia Syringes

Suitable For Use with All 1.8ml Dental Catridges






MS7-101



MS7-102






MS7-103



MS7-104


Mechanical Aspirating Local Anesthesia Syringes

Suitable For Use with All 1.8ml Dental Catridges






MS8-101



MS8-102





MS8-103



MS8-104

Intraligamental Syringes



MS9-101



MS9-102



MS9-103

Made Of Stainless Steel Color Coated. Suitable For Use with All 1.8mlDental Cartridges. Inject 0.05ml with Every Click. Exact Dosage Every Time.Supplied with Angled Tubes.

Gun Style Intraligamental Syringes


Use Mr shah Premium Quality Stainless Steel Syringe togivethe anesthesia to patients without any problem.



MS9-104 STRAIGHT
MS9-105 ANGLED
Made of Stainless Steel
Mr Shah Ject
Metric Thread
Inject 0.06ml with every click
Exact dosage every time
Suitable for use with all 1.8ml dental Cartridges
Supplied with Straight and angled tubes

URL: www.mrshah.co
TEL:+49 174 8102485 | Email: info@mrshah.co

Handles for Mouth mirrors



MS10-101 Mirror Handle Stem Round (6mm) Length 13cm

MS10-102 Mirror Handle Stem octagonal (6mm) Length 13cm

MS10-103 Mirror Handle Slim Shape Pattern (6mm) Length 15cm

MS10-104 Hollow Mirror Handle Stem Extra Light (7.5mm) Length 13cm

MS10-105 Mirror Handle Stem with Rule 35mm (6mm) Length 13cm

MS10-106 Mirror Handle Stem with Rule 35mm (6mm) Length 13cm **(Cone Socket)**

MS10-107 Hollow Mirror Handle Stem Extra Light (7.5mm) Length 13cm **(Cone Socket)**

MS10-108 Hollow Mirror Handle Slim Shape Pattern (9mm) Length 13cm **(Cone Socket)**

MS10-109 Hollow Mirror Handle Stem Extra Light (9mm) Length 13cm **(Cone Socket)**

MS10-110 Hollow Mirror Handle blue Titanium Coated (9mm) Length 13cm

MS10-111 Hollow Mirror Handle Light weight(9.5mm) Length 13cm **(Cone Socket)**

Dental Mouth Mirrors



Mouth Mirror plain



MS11-101 Fig. 3 ∅20mm

MS11-102 Fig. 4 ∅22mm

MS11-103 Fig. 5 ∅24mm



Mouth Mirror Cone Socket Front Surface



MS11-104 Fig. 3 ∅20mm

MS11-105 Fig. 4 ∅22mm

MS11-106 Fig. 5 ∅24mm

Double Sided Mouth Mirror

Front Surface



MS11-107 #4 (22mm)
MS11-108 #5 (24mm)

MS11-109 #4 (22mm)
MS11-110

Matrix Retrainer



MS12-101 Ivory 63mm



MS12-102 Tofflemire-Universal



MS12-103 Tofflemire-Junior



MS12-104 Ivory8 5mm

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

Dental Tweezers

- Tweezers are used particularly to take and place cotton pellets, wedges and articulation
- Tweezers are important Dental Lab Equipment, made of non-magnetic stainless steel
- They are designed particularly to hold the tiniest article, provide better grip on the ja paper
- from the to the mouth and back.





MS13-101 London College 15cm
MS13-102 Meriam 16cm
MS13-103 London College 15cm
MS13-104 Meriam 16cm
MS13-105 14cm 1x2 Tissue Forcep
MS13-106 College with Lock 15cm cvd
MS13-107 Hamingway 15cm
MS13-108 Adson Brown 12cm
MS13-109 Adson 12cm
MS13-10 Adson 12cm 1x2



MS13-111 Bonding Bracket Holder 13.5cm
MS13-112 Bonding Bracket Holder 13.5cm
MS13-113 Bonding Bracket Holder 13.5cm
MS13-114 Miler 15cm Str.
MS13-115 Miler 15cm cvd
MS13-116 Bonding Bracket Holder 13.5cm
MS13-117 Bonding Bracket Holder 13.5cm
MS13-118 Bonding Bracket Holder 13.5cm
MS13-119 Bonding Bracket Holder 13.5cm
MS13-120 Corn Suture Pliers, 15cm

Adson Tweezers (Tungsten Carbide)




MSTC1-101 Adson 12cm 1x2
MSTC1-102 Adson 12cm
MSTC1-103 Adson, 12cm

Micro Tweezer



MSTC-104 Micro Tweezer 17cm

Lab Tweezers



MSLT-101 Forceps 14cm Curved, Ergonomic Paddle Handle, 1.5mm tip, smooth.
MSLT-102 Forceps 14cm Straight, Ergonomic Paddle Handle, 1.5mm tip, smooth.
MSLT-103 Forceps 14cm Curved, Ergonomic Paddle Handle, 1.5mm tip, serrated.
MSLT-104 Forceps 14cm Straight, Ergonomic Paddle Handle, 1.5mm tip, serrated.
MSLT-105 Forceps 14cm Straight, Ergonomic Paddle Handle, Needle tip.

URL: www.mrshah.co
TEL:+49 174 8102485 | Email: info@mrshah.co

Functions and Features:

- Scissors designed for suture removal •
- Blades - one has a curved edge and one has a notched edge
- During suture removal, suture is scooped up by blade with the3.0 mm notched end before cutting to prevent tissue trauma

General Scissors



MS14-101
Iris, 11.5cm
Straight

MS14-102
Iris
11.5cm
Curved



MS14-103
Neumann
12cm
Double Curved



MS14-104
Stevens
13cm,Straight

MS14-105
Stevens
13cm,Curved



MS14-106
Locklin 16cm
Angular



MS14-107
Spencer
13cm



MS14-108
Operating Scissor
13cm Straight

MS14-109
CVR

Castroviejo Scissors

Extra Light Weight





Curved
MSTT1-103
15cm
MSTT1-104
17cm

MSTT1-101 - 15cm
MSTT1-102 - 17cm
Castroviejo
Scissor, Straight

Made By Titanium

Mr Shah Flex-Spring Scissors



Curved
MS14-112
15cm
MS14-113
17cm



Curved
MS14-115



MS14-110-15cm
MS14-111-17cm
Castroviejo
Scissors, Straight

MS14-114
Noyes
11cm
Straight

MS14-116
Noyes
11cm
Angular

Super Cut Scissors



MS14-118
Curved



MS14-120
Curved



MS14-122
Curved

MS14-117
Iris
11.5cm Straight
Saw edge

MS14-119
Kelly
16cm Straight
Saw edge

MS14-121
Beebee 11cm
Saw edge, Straight

Scissors Tungston Carbide


Tungsten Carbide scissors feature cutting edges of permanently welded Tungston carbide inside. These Scissors will retain a fine sharp cutting edge for long time and will outlast any standerd stainless steel dental scissors. Both handles are gold-plated for easy identification.





Curved
MSTC2-103
15cm
MSTC2-104
17cm

MSTC2-101-15cm
MSTC2-102-17cm
Castroviejo
Scissors, Straight



MSTC2-106
Iris
11.5cm
Curved

MSTC2-105
Iris
11.5cm
Straight





MSTC2-108
Goldman-Fox
13cm
Curved

MSTC2-107
Goldman-Fox
13cm
Straight



MSTC2-110
Kelly
16cm
Curved

MSTC2-109
Kelly
16cm
Straight





MSTC2-112
Beebee
11cm
Curved


MSTC2-111
Beebee 11cm
Saw edge, Straight



MS14-123
Universal
12cm
Saw edge

MR SHAH Needle holders feature tungsten carbide inserts with a matched diamond-cut pyramidcal profile for enhanced longevity, a secure, positive suture needle grasp and reduced risk of cutting the suture material. With a precision designed hings, the smooth surface reduces the risk of catching and or cutting the suture material during the knot tying process, enhancing clinical efficiencies. For optimal grip and control, **MR SHAH** offers a variety of design options, including standard, finger and palm controlled closure.

Tungsten Carbide Needle Holders



MSTC3-101 Baumgarte 13cm

Mayo-Hegar **MSTC3-102** 130mm **MSTC3-103** 160mm

Olsen-Hegar Needle Holder With Scissors **MSTC3-106** 12cm

Castroviejo **MSTC3-109** STR 14cm **MSTC3-110** STR 16cm **MSTC3-111** STR 18cm

Castroviejo Fine Pointed **MSTC3-115** STR 14cm **MSTC3-116** STR 16cm **MSTC3-117** STR 18cm

MSTC3-121 Mathieu Needle Holders 14cm.

Needle Holders




MS15-101 Halsey 13cm

MS15-102 Crile-Wood 15cm

MS15-103 Mathieu Needle Holders 14cm

Titanium Made Instruments

Lightweight With A High Tensile Strength




MSTT2-101 Halsey 13cm

Castroviejo Fine Pointed **MSTT2-102** STR 14cm **MSTT2-103** STR 16cm **MSTT2-104** STR 18cm

MSTT2-108 Mosquito Forceps 12.5cm,Curved

Forceps



MS16-101 Hartmann (Baby Mosquito) 9cm Straight

MS16-103 Kelly-Rankin 16cm straight

Halstead Mosquito **MS16-105** 12.5cm straight

Halstead Mosquito 1x2 **MS16-107** 12.5cm straight

Loma **MS16-109** 10cm **MS16-110** 13cm

MS16-111 Allis Tissue Forceps 15cm

Foerster Dressing Forceps **MS16-112** 20cm


A1-Double Ended


A2-Double Ended


A3-Double Ended


A4-Double Ended


A5-Double Ended


A11-Double Ended


A6-Double Ended


A7-Double Ended


A8-Single Ended


A9-Single Ended


A10-Single Ended


A12-Double Ended


A13-Single Ended

HOLLOW HANDLES (A)

TITANIUM COATED HANDLES (B)


B1-Double Ended


B2-Double Ended


B3-Double Ended


B4-Double Ended


B5-Double Ended


B6-Double Ended



B7-Single Ended



B8-Single Ended



SOLID HANDLES (D)


D1-Double Ended


D2-Double Ended


D3-Double Ended


D4-Double Ended


D5-Single Ended


D6-Double Ended

SILICON HANDLES (E)


E1-Double Ended


E2-Double Ended



E3-Double Ended


E4-Double Ended



E5-Double Ended


E6-Double Ended


Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)



Silicon Handls (E)




Explorers Double Ended



MEX-101
Explorers Double Ended 2A


Probes Double Ended



MEX-110
Probe Double Ended 12


EXPROS

A double-ended instrument that combines two critical instruments required for every effective diagnostic set-up with an explorer on one end and a probe on the other, an Expros is excellent for dental examinations and charting.



Excavators


SPOONS: For Removal of Carious Dentin



MEC-101-E1

Pluggers / Condensers


Used to Pack Filling Materials into Cavity Preparations




MPL-110- 1/2,ANDREW

CARVERS

Used to Crave anatomical features and trim excess materials.




Hollow Handls (A)




Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Burnishers


Designed to condense, smooth, carve and polish amalgam.
Gold Titanium Coated Points.



MBU-101-18

PK-Thomas

Used for waxing procedures and techniques.



MPKT-101 **MPKT-102** **MPKT-103** **MPKT-103R** **MPKT-104** **MPKT-105**

Carving Knives

Used for trimming excess filling material, flash and overhangs.



MCK-113 / 20 Esthetic
For anterior teeth. Sharp, thin blade allows access to all surfaces. The offset angle provides universal adaptability



MCK-114 / 21 Esthetic
For posterior teeth. Thin, sharp offset angle provides access to surfaces where linear finishing strips would not be effective



MCK-115 / 6 Tanner
For posterior teeth. Sharp offset angle provides access to many surfaces



MCK-116 / 12 Novatech Cement Remover
Combines a sickle-shaped scaler with a flat blade for removal of excess resin. provides access to many surfaces cement. or porcelain flasn
The narrow chisel removes excess interproximal material with a push stroke



MCP-116 / BN1
Thin blade and rounded contour facilitates use in both thick and thin tissues without catching or dropping cord. Bilateral notch allows placement in limited access areas.



MCP-117 / 1 Yardley
Round non-serrated working end.



Serated **MCP-118** Non-Serated
CSI-1



MCP-119
7 Guyer



Serated **MCP-120** Non-Serated
113

Gingival Retractors

Protects tissue during cavity preparations and placement of composite restorations




MGR-120 / 1 Meinershagen / For mandibular inscisors



MGR-121 / 2 Meinershagen
For maxillary and mandibular premolars and canines



MGR-122 / 3 Meinershagen
For maxillary central incisors and wide canines



MGR-123 / 4 Meinershagen
For all molars.



MGR-124
10 Goldman-Fox


Composite Filling Instruments

Also Available in Titanium Coated Colored Points.

Vertical and horizontal combination.




MCF-101



MCF-102

MCF-103




MCF-104

MCF-105



MCF-106

MCF-107



MCF-108

MCF-109




MCF-110
Woodson




MCF-111
W1

Vertical in and offset angled.



MCF-112


MCF-113



MCF-114

MCF-115


Hollow Handls (A)



Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Composite Filling Instruments



MCG-101 /#1 Small Universal



MCG-102
#2 Medium Universal



MCG-103
#4 Small Universal



MCG-104
#3 Large Universal



MCG-105
#5 Medium Universal



MCG-106-Boghosian
Combination of Thin, Knife-Shaped Blade



MCG-107-8A
Use for packing gingival retraction cord.



MCG-108-Boghosian
Used for measuring composite layers.



MCG-109-Contouring
Used for shaping of inclines, planes
for Anterior and posterior restorations.



MCG-110
Quétin



MCG-111-8A
#2 Round Arrow

Composite Filling Instruments Set



MCB-125



MCB-126



MCBS1-01
Set of 6 in Cassette



MCB-127



MCB-128




MCB-129



MCB-130

Amalgam Carrier



MS15-110
Mini/Regular 90°


Amalgam Well



MSB-10D



MS15-111
Jumbo/Regular 30°
20 cm



MS15-112
Large/Regular 45°
20cm



MS15-113
mm,2,00 ø
15cm



MS15-114
mm,2,00 ø
15cm



MS15-116
mm,2,00 ø
15cm



MS15-116

Sickle Scalers

Designed for removal of supragingival calculus.

All Scalers & Cutters are Available in

J2-420 ECO	USA 440c Premium
Best Quality Medical Stainless Steel. High Chromium Content Good Strength Easy To Sharpe Economical Cost Available in Different Sizes	440c High-end Stainless Steel High Chromium Content High Carbon Martensitic High Strength Better Edge Retention Titanium Coated Points



MSC-101
Contra angle design offers greater access to distal surfaces



MSC-102 / H6/H7
Super-thin, contraangle design is for interproximal areas. Shank angulation provides accessibility to all anterior and premolar surfaces



MSC-103
H5/U15
Hygienist-Towner



MSC-104 / H5/33 - Hygienist-Jacquette
Super-thin sickle scaler combined with a straight sickle blade. For scaling interproximal and cervical areas of anterior teeth



MSC-105
U15/30
Towner-Jacquette



MSC-106
U15/33
Towner-Jacquette



MSC-107
204



MSC-108
For posterior supra-gingival interproximal scaling



MSC-109 - 0/00 Morse
Miniature blade for scaling crowded or overlapping anterior teeth excellent for stain removal



MSC-110/204SD
Contra angle design offers greater access to distal surfaces



MSC-111
204S



MSC-112
Nabraska 5/6



MSC-113
NV1/H5



MSC-114
Nevi-1 / H5
Anterior



MSC-115
Nevi-2 / Posterior
A super thin posterior curved sickle scaler

Columbia Curettes



MCC-106
13/14 Columbia
(Regular)

Universal Application

MCC-107
13/14 Columbia
(Rigid)




MCC-108
2R/2L Columbia, Anterior application



MCC-109
4R/4L Columbia, Posterior application


Hollow Handls (A)



Titanium Coated Handls (B)




Solid Handls (D)



Silicon Handls (E)




Hollow Handls (A)




Titanium Coated Handls (B)



Solid Handls (D)



Silicon Handls (E)



Langer Curettes



MLC-101
1/2 Langer
(Reg-shown)

MLC-102
1/2 Langer
(Rigid)

Langer 1/2 for
mandibular posterior

MLC-103
1/2 Langer
After Five

MLC-104
1/2 Langer
Mini Five



MLC-105
3/4 Langer
(Reg-shown)

MLC-106
3/4 Langer
Rigid

Langer 3/4 for
maxillary posterior

MLC-107
3/4 Langer
After Five

MLC-108
3/4 Langer
Mini Five



MLC-109
5/6 Langer
(Reg-shown)

MLC-110
5/6 Langer
Rigid

Langer 5/6 for
maxillary and
Mandibular anterior

MLC-111
5/6 Langer
After Five

MLC-112
5/6 Langer
MiniFive



MLC-113
17/18 Langer
(Reg-shown)

Langer 17/18 for maxillary
and Mandibular posterior

MLC-114
17/18 Langer
Rigid

Barnhart



MBC-105
1/2 Barnhart
(Reg-shown)

MBC-106
1/2 Barnhart
(Rigid)



MBC-107
5/6 Barnhart
(Reg-shown)

MBC-108
5/6 Barnhart
(Rigid)

Periodontal Files



MPF-150
9/10 Schluger
Curved



MPF-151
11/12 Buck
Straight



MPF-152
1S/2S Sugarman
Mesial/Distal



MPF-153
3S/4S Sugarman
Buccal/Lingual
For removal of granulonatous tissue in
supra-bony and infrabony pockets.

- Area-specific to allow for deep scaling, root planing and periodontal debridement
- Offset blade provides a perfect working angulation for the tooth surface
- Offset blade angulation means only the lower cutting edge of each blade is used

Gracey Standard Curettes





MGC-110 1/2 (Reg-shown)	**MGC-110-1** 1/2 (Rigid)	Slight contra-angle for anterior incisors and canines	**MGC-110-2** 1/2 After Five	**MGC-110-3** 1/2 Mini Five



MGC-111 3/4 (Reg-shown)	**MGC-111-1** 3/4 (Rigid)	Short contra-angle for anterior Incisors and canines	**MGC-111-2** 3/4 After Five	**MGC-111-3** 3/4 Mini Five



MGC-112 5/6 (Reg-shown)	**MGC-112-1** 5/6 (Rigid)	Medium contra-angle for anteriors and premolars	**MGC-112-2** 5/6 After Five	**MGC-112-3** 5/6 Mini Five



MGC-113 7/8 (Reg-shown)	**MGC-113-1** 7/8 (Rigid)	Medium contra-angle for premolars and molas, facial and lingual surfaces	**MGC-113-2** 7/8 After Five	**MGC-113-3** 7/8 Mini Five



MGC-114 9/10 (Reg-shown)	Long contra-angle for molar buccal and lingual surfaces and a more oronounced angle than 7/8 for less accessible root surfaces.	**MGC-114-1** 9/10 (Rigid)	



MGC-115 11/12 (Reg-shown)	**MGC-115-1** 11/12 (Rigid)	Angulated to reach mesial surfaces of posterior teeth	**MGC-115-2** **11/12** After Five	**MGC-115-3** **11/12** Mini Five



MGC-116 13/14 (Reg-shown)	**MGC-116-1** 13/14 (Rigid)	Angulated to reach distal surfaces of posterior teeth	**MGC-116-2** **13/14** After Five	**MGC-116-3** **13/14** Mini Five




MGC-117 15/16 (Reg-shown)	**MGC-117-1** 15/16 (Rigid)	Same shank bends as 13/14, but blade positioned to reach mesial posterior surfaces	**MGC-117-2** **15/16** After Five	**MGC-117-3** **15/16** Mini Five



MGC-118 17/18 (Reg-shown)	Accentuated angles improve access to distal posterior surfaces. The 17/18 also has smaller blade and silany longer terminal shank.	**MGC-118-1** **17/18** (Rigid)


Hollow Handls (A)




Titanium Coated Handls (B)




Solid Handls (D)




Silicon Handls (E)




Hollow Handls (A)



Titanium Coated Handls (B)



Solid Handls (D)




Silicon Handls (E)



Periosteals

For reflecting and retracting the mucoperiosteum after incision of the ginbival tissue.



MPS-101
9 Molt

MPS-102
P9/PA

MPS-103
1 Hopkins

MPS-104
Demirel

MPS-105
3 Mead

MPS-106
8D

MPS-107
9H Howard

MPS-108
9A Allen

MPS-109
24G

MPS-110
15 Freer Curved

MPS-111
Simion

MPS-112
152 Kramer Nevins Blunt

MPS-113
PR-3 Prichard

MPS-114
3/9A

MPS-115
#3M Curved

MPS-116
1 Lahanca

MPS-117
2 Labanca Modified

MPS-118
Raspator, Prichard

MPS-119
Mixing Spatula

Goldman-Fox



MGF-154-8 Goldman-Fox
Spear Point - For cutting interdentally after initial incision. Severs buccolingually. Angulated blade for easier handling.



MGF-155-9 Goldman-Fox
Combination - Spear-pointed blade used in mandibular region from posterior direction. Modified sickle blade (elongated) used in interdental area for debridement.



MGF-156-11 Goldman-Fox
Spear Point - Similar to #8, except blades are offset for use in posterior region.



MGF-157
3/4 Buck



MGF-158
5/6 Buck

Retractors & Mouth Gags



SM17-101-Infant, Small
SM17-102-Childern, Medium
SM17-103-Adults, Large



Weitlaner
SM17-104-10.5cm
SM17-105-13.5cm



SM17-106
Stainless Steel
Double Ended Cheek Retractors



SM17-107
Stainless Steel
Single Ended Cheek Retractors



Cheek Retractors,
SM17-108-14cm
SM17-108-18cm



Minnesota
SM17-110-Small
SM17-111-Large



SM17-112
Retractors
17cm
Hooks 22 mm

SM17-113
Retractors,Cheek
retractors combination
21 cm,2



SM17-114
140 mm

SM17-115
150 mm



SM17-116
Bishop



SM17-117
Austin Fig.1




SM17-118
Austin Fig.2



Surgical Retractor
SM17-119-7x25mm
SM17-120-10x42mm
SM17-121-12x55mm

SM17-122
Implant Retractor
23cm



MSAS-01
6 Way Adjustable Angled Scalpel Handle

Scalpel Handles



SM18-101 No. 5E
SM18-102 No. 5A
SM18-103 3K Mini
SM18-104 Fig.3
SM18-105 Fig.4
SM18-106 Fig.3G

Bone Files



SM19-101 Miller Cross Cut
SM19-102 Miller 2X Straight Cut
SM19-103 Wahi 1


Cement Spatulas & Wax Knives



SM20-101 Gritman
SM20-102 3
SM20-103 1 Medium Cement
SM20-104 Small **SM20-105** Large Wooden Handle Wax Knives,
SM20-106 Wooden Handle Fig. 8 20 cm curved

TEL:+49 174 8102485 | Email: info@mrshah.co URL: www.mrshah.co

Rubber Dam Forceps



MSS21-01
01-Anisworth
02-Brewer
03-Clamps Stand
04-Rubber Dam Frame
05-Rubber Dam Clamp
Fig,212,211,8,8A
14,0,00.2,2A,7



MSS21-02
01-Anisworth
02-Palmer
03-Clamps Stand
04-Rubber Dam Frame
05-Rubber Dam Guide
06-Rubber Dam Clamp
Fig,212,211,8,8A
14.0.00.2.2A.7

RUBBER DAM FRAME
MS21-107
Young
Adult / Child.





MS21-101
Ivory
160mm



MS21-102
Anisworth

MS21-103
Brewer



MS21-104
Ivory light weight



MS21-105
Palmer



MS21-106
Seaparating Elastic Plier

Rubber Dam Clamps


 **MS21-108** Fig.15	 **MS21-109** Fig.18	 **MS21-110** Fig.24	 **MS21-111** Fig.25	 **MS21-112** Fig.26	 **MS21-113** Fig.28	 **MS21-114** Fig.30	 **MS21-115** Fig.31	 **MS21-116** Fig.51	 **MS21-117** Fig.138
 **MS21-118** Fig.139	 **MS21-119** Fig.200	 **MS21-120** Fig.201	 **MS21-121** Fig.202	 **MS21-122** Fig.203	 **MS21-123** Fig.204	 **MS21-124** Fig.205	 **MS21-125** Fig.3	 **MS21-126** Fig.W3	 **MS21-127** Fig.4
 **MS21-128** Fig.56A	 **MS21-12+** Fig.7	 **MS21-130** Fig.7A	 **MS21-131** Fig.8	 **MS21-132** Fig.8A	 **MS21-133** Fig.W8A	 **MS21-134** Fig.10	 **MS21-135** Fig.11	 **MS21-136** Fig.12A	 **MS21-137** Fig.13A
 **MS21-138** Fig.14	**MS21-139** Fig.14A	**MS21-140** Fig.14AN	**MS21-141** Fig.22	**MS21-142** Fig.27	**MS21-143** Fig.207	**MS21-144** Fig.29	**MS21-145** Fig.206	**MS21-146** Fig.208	**MS21-147** Fig.209
MS21-148 Fig.1	**MS21-149** Fig.1A	**MS21-150** Fig.2	**MS21-151** Fig.2A	**MS21-152** Fig.W2A	**MS21-153** Fig.0	**MS21-154** Fig.00	**MS21-145** Fig.210	**MS21-156** Fig.211	**MS21-157** Fig.212
MS21-158 Fig.6	**MS21-159** Fig.9	**MS21-160** Fig.N06	**MS21-161** Fig.U67	**MS21-162** Fig.214	**MS21-163** Fig.P1	**MS21-164** Fig.P2	**MSS21-03** Rubber Dam Clamp 9 pcs Kit Fig.7,14A,2A,8A,4 13A,12A ,9,2		



MS22-101
Dental
crown remover
forceps



MS22-102
Crown Spreaders



MS22-103
WYNMAN
Crown Gripper,CARB-N-SERT



MS22-104
Key Crown Remover
Set of 3 pcs



MS22-105
Universal
Crown Remover
Hollow Handle



MS22-106
Automatic
Crown Remover
with 2 Tips

Bone Gauges and Measuring Instruments



MSS22-01
Crown Remover Set
Full Automatic



MSS23-101
Boley Gauge



MS23-102



MS23-103

Bracket Positioning Boons Gauge, (Star Gauge)

MS23-107

- It Is The Most Frequently Used Bracket Positioning Aid.
- It is Made Of Stainless Steel.
- It Accurately Measures Height Of Bracket Placement From The Incisal Edge At 3.5, 4, 4.5 And 5mm.
- Flat Surface Rests On Incisal/occlusal Surface Of Tooth. Simplife. Seating The Bracket On The Teeth With Exactly The Right Distance Between Bracket Dot And Incisal Edge With Pencil Lid Or Metal Tips.





MS23-104
Measuring Loop
Straight



MS23-105
Mearuring Loop
Curved

MS23-106
Lwanson Crown
Guage




MS22-107
Napkin Holder
Silicon

MS22-108
Napkin Holder

Bone Pluggers



Sinus Lift Instruments



Bone Chisels




MS26-102 8mm

MS26-104 8mm

MS26-101 Kramer-Nevins Bone chisel 5mm

MS26-103 Ochsenbien 5mm

Tissue Punch Straight & Curved




MS27-111 2.0mm
MS27-112 4.0mm
MS27-113 5.0mm
MS27-114 6.0mm
MS27-115 7.0mm
MS27-116 8.0mm

MS27-118 2.0mm
MS27-119 4.0mm
MS27-120 5.0mm
MS27-121 6.0mm
MS27-122 7.0mm
MS27-123 8.0mm

MS27-110 Tissue Punch Stright

MS27-117 Tissue Punch Curved

Osteotome Concaves & Concaves



MSS28-01 Set Of 5 Pieces Concaves Straight
MSS28-02 Set Of 5 Pieces Convexes Straight

MSS28-03 Of 5 Pieces Concaves Curved
MSS28-04 Of 5 Pieces Convexes Curved



MS29-101
Bone Ridge Mapping Forceps



MS29-102
Measuring Caliper



MS29-103
Splinter Forceps Straight

MS29-104
Splinter Forceps Curved



MS29-105
Bone Mill Plier



MS29-106
Molt 18.5cm Auto Clavable Plastic Tips.

Bone Scraper



MS30-101 Curved
MS30-102 Straight

Star Torque Wrench



MS31-101
Star Torque Wrench 6.35mm Universal



MS29-108
Bone Grafting Syringe
17 cm straight / Bend
I.D 3.5mm, O.D 4.5

MS29-109
Bone Grafting Syringe
18 cm straight / Bend
I.D 4.5mm, O.D 5.5

MS29-110
Bone Grafting Syringe
18 cm straight / Bend
I.D 6.0mm, O.D 7.0

MS29-107
Bone Grafting Syringe
17 cm, Straight / Bend
I.D 2.5mm, O.D 3.5

MS29-111
Suction Tub



MS29-113
Ø 3mm
Length 180mm
MS29-114
Ø 4mm
Length 180mm
MS29-115
Ø 5mm
Length 180mm

MS29-112
Fig. 6FR Frazier
15 cm-ø2mm



MS29-116
Debakey Atraumatic Tissue Forceps Str Serrated Teeth 16cm

MS29-117
Debakey Atraumatic Tissue Forceps Cvd Serrated Teeth 16cm

Bone Expander Kit



MSBE-01
Bone Expander Kit
Expander:2.6mm, 3.0mm
3.4mm, 3.8mm, 4.3mm
Rachet Driver
Parallel Pin
Pilot Drill
Drill 1.4mm
Drill Driver 4mm
Hand Driver 4mm
Rachet 4.0mm



MS29-118
Bone Crusher Well

Periotomes



MS32-101 / Pt5 Titanium TG CVD S/E
MS32-102 / Pt4 Titanium STR S/E
MS32-103 / Pt3 Titanium Angle S/E
MS32-104 / Pt1 Titanium A/V Regular
MS32-105 / Pt2 Titanium A/V Regular

Micro Serrated



MS32-120
Mini Drill Set of 8 Pcs.
in Stainless Steel Stand



MS32-121
Standard Drill Set of 8 Pcs.
in Stainless Steel Stand



MS32-122
Trephine Bur Set of 8 Pcs.
in Stainless Steel Stand



MS32-123
Bone Compression
Prosthodontic Dental Kit
for Implant



MS32-124
Tissue Punch KiT
inStainless Steel Stand



MS32-127
bur holder
for 12 Pcs

MS32-128
bur holder
for 06 Pcs



MS32-125
Trimmer Kit in
Stainless Steel Cassette



MS32-126
Mini Expander Kit in
Stainless Steel Cassette

Soft Tissue Brushing Instruments

MSTB-Set of 5 Pcs



MSTB-1
Large Tip



MSTB-2
Medium Tip



MSTB-3
Regular Small Tip



MSTB-4
Angulated Right Tip (Posterior)



MSTB-5
Angulated Left Tip (Posterior)



MS32-129
Conical Universal Kit
(Kit of 50 PCS)
Stainless Steel Cassette

URL: www.mrshah.co
TEL:+49 174 8102485 | Email: info@mrshah.co

Distal End Cutters



Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)

MS33-101
Distal End
Cutter Universal Cut
long handle 14cm



MS33-102
Distal End Cutter
Universal Cut
long Neck 12cm

Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)



MS33-103
Distal End Cutter
Universal Cut
long handles 14cm

Max: .021"(.54mm)
x.025"(.64mm)



MS33-104
Flush Cutter with
Hold 12cm

Max: .021"(.54mm)
x.028"(.71mm)
Min:.012"(.31mm)



MS33-105
Distal End Cutter
Universal Cut
long handles 14.5cm

Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)



MS33-106
Distal End Cutter
Universal Cut
Box Joint 12cm

Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)



MS33-107
Distal End Cutter
Universal Cut
long neck 12cm

Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)



MS33-108
Distal End Cutter
Universal Cut heavy
pattern,Box Joint long
Handles 14cm

Max: .021"(.54mm)
x.025"(.64mm)
Min:.012"(.31mm)

Ligature & Pin Cutters



Max: .021"(.31mm)
Soft Wires

MS33-109
Ligature & Pin Cutters
straight 12.5cm



MS33-110
Mini Ligature & Pin
Cutters 15 degree
angled, Slim Beak
12.5cm

Max: .012"(.31mm)
Soft Wires



MS33-111
Ligature & Pin
Cutters wide beak
15 degree angled
12.5cm

Max: .012"(.31mm)
Soft Wires



MS33-112
Ligature & Pin
Cutters straight
12cm

Max: .012"(.31mm)
Soft Wires



MS33-113
Micro Ligature & Pin
Cutters 15 degree
angled 11.5cm

Max: .012"(.31mm)
Soft Wires



MS33-114
Ligature Cutters
45 degree angled
wide Handles
Very Fine Tips.
11.5cm

Max: .012"(.31mm)
Min: .008"(.20mm)
Soft Wires



MS33-115
Ligature & Pin
Cutters, Box
Joint 1/3 Handles
Gold Platted 12.5cm

Max: .021"(.31mm)
Soft Wires



MS33-116
Micro MiniLigature
Cutters 15 degree
angle 11cm

Max: .021"(.31mm)
Soft Wires

Pliers



MS33-117
Bracket Removing
Plier 12.5cm angled



MS33-118
Bracket Removing
Plier 12.5cm
straight



MS33-119
Band Removing
Plier long neck
12.5cm



MS33-120
Slim Weingart
Plier 13cm



MS33-121
Weingart
Plier 12.5cm



MS33-122
Weingart
Plier 14cm



MS33-123
Nance Pliers
12.5cm



MS33-124
Loop Forming
Pliers 13cm



MS33-125
Three Jaw Pliers
12.5cm

Max: .028"(.71mm)
Wires



MS33-126
Adrere Pliers
13cm

Max: .030"(.76mm)
Wires



MS33-127
Bird Beak
Pliers 12cm

Max: .030"(.76mm)
Wires



MS33-128
Lab Bird Beak
Pliers with serrated
tips 14cm

Max: .036"(.92mm)
Hard Wires



MS33-129
Light Wire Bending
Pliers long 12.5cmted

Max: .020"(.51mm)
Wires



MS33-130
Slim Bird Beak
Pliers Tapered
Beak 12.5cm

Max: .022"x0.25"
(.56mmx.64mm)
Wires



MS33-131
Rectangular Arch
Forming Pliers
12cm

Max: .022"x0.25"
(.51mmx.64mm)
Wires



MS33-132
Slim Arch
Forming Pliers
12cm

Max: .022"x0.25"
(.51mmx.64mm)
Wires



MS33-133
Tweed Arch
Forming Pliers
12cm

Max: .022"x0.25"
(.51mmx.64mm)
Wires



MS33-134
Hollow chop
Arch Forming
Pliers 12cm



MS33-135
Hollow chop
Arch Forming
Pliers with 3
Grooves 12cm



MS33-136
Tweed Pliers
12.5cm



MS33-137
Adhesive Removing
Pliers



MS33-138
How Utility
Pliers angled
12.5cm



MS33-139
How Utility
Pliers Straight
12.5cm



MS33-140A
Tweed Pliers
12.5cm

Max: .022"x0.28"
(.56mmx.71mm)
Wires



MS33-140B
Loop Forming
Pliers with cutters.
Tweed Type

Max: .025"
Rectangular Hard wires



MS33-141
Step Forming Pliers
A-0.25 / **B**-0.50
C-0.75 / **D**- 1mm

Max: .022"
(.56mm) Wires



MS33-142
Ligature Forming
Pliers



MS33-143
Niti Cinch Back
Pliers, Slim Beak

Max: .021"x0.25"
(.54mmx.64mm)
Wires



MS33-144
Hook Cimping
Pliers



MS33-145
Loop Forming
Pliers,long Beak

Max: .028"(.70mm)
Hard Wires



MS33-146
Optical Pliers



MS33-147
Young Loop
Bending Pliers

Max: .028"(.70mm)
Hard Wires



MS33-148
Crown Crimping
Pliers



MS33-149
Johson Contouring
Pliers



MS33-150
Niti Cinch Back
Pliers, 13.5cm
Short Handle
Max: .021"x0.25"
(.54mm x .64mm)
Wires

Clear Brace Aligner Pliers



MS34-201
Hole Punch
Cutter



MS34-202
Tear Drop
Cutter



MS34-203
Horizontal
Dimple



MS34-204
Mini Dot
Dimple




MS34-205
Vertical Dimple



MS34-206
Dimple Remover

Individual Torquing Pliers Set



MS33-155
Individual Torquing Pliers set
To make 6mm long torque 13cm

Can bend wire upto. .022" x.025"

Hard Wire Cutters



MS36-101
Hard wire Cutter



MS36-102
Hard wire Cutter



MS36-103
Hard wire Cutter



MS36-104
Hard wire Cutter


WIRE CUTTER TC GOLD

Tungsten Carbide (Tip Insert)



MS36-01-16 cm
NS36-02-18 cm
wire Cutter
TC Gold
Double Action

Hand Instruments



Ortho light scaler double ended — **MS36-105**
Twiel-on elastic module instrument — **MS36-106**
Ortho Light double ended ligature director — **MS36-107**
Ortho Light double ended ligature scaler/director tivate coil spring are remove modules with one instruments — **MS36-108**
Schure Band Seater — **MS36-109**
Band seater auto cleavable — **MS36-110**
Mershon Band Seater — **MS36-111**

Uses of Articulator

Primary purpose:

1) To hold opposing casts in a predetermined fixed relationship
2) To open and close. 3) To produce border and inter- border diagnostic sliding motions of the teeth similar to these in the mouth.







MSAR-1001
Technical Datas:
Intercondylar distance 120 mm. / condylar inclination fixed at 30°
Bennet angle fixed at 7° / Incisal table inclination 10°
Distance between mounting plates 85 mm. / W eight 400 gr.
Opening 180°

MSAR-1002



Black

MSAR-1003



Blue

MSAR-1004



Yellow

TEL:+49 174 8102485 | Email: info@mrshah.co
URL: www.mrshah.co

Impression Trays Dentulous Perforated

MS37-101



SIZE: S-M-L-XL

Impression Trays Edentulous Perforated

MS37-103



Impression Trays Dentulous Perforated

MS37-102



No Welding
Easy to Clean
Adjustable Front Wall

SIZE: S-M-L-XL

Impression Trays Edentulous Non-Perforated

MS37-104



SIZE: S-M-L-XL

Impression Trays Dentulous Perforated

MS37-110



No Welding
Easy to Clean
Adjustable Front Wall

SIZE: S-M-L

Impression Trays Edentulous Perforated

MS37-111



No Welding
Easy to Clean

SIZE: S-M-L

Impression Trays Dentulous Non-Perforated

MS37-112



SIZE: S-M-L

Mrshah Impression Trays Size Chart

Uppers / superiori Sizes / Misure (mm)

	U1 XXS	U2 XS	U3 S	U4 M	U5 L	U6 XL
	56	65	68	71	73	78
	44	62	54	58	59	63

Lower Impression Trays Sizes Chart

	L1 XXS	L2 XS	L3 S	L4 M	L5 L	L6 XL
	56	66	70	73	76	78
	44	50	52	56	58	60

All Mr Shah Impression Trays are Made by Surgical Grade Stainless Steel

Impression Tray Half
Perforated



MSH-101
U1P / 12,5cm
For upper left jaw
For lower right jaw

MSH-102
U1P / 12,5cm
For upper right jaw
For lower left jaw



MSH-103
U1P / 12,5cm
Adjustable Perforated



MSH-104
U1P / 12,5cm
Adjustable Solid



MS38-101
Cotton Holder



MS38-102
Spirit Lamp



MS38-103
Endodontic Box



Kidney Trays

With lid available	Without lid available
MS38-104 6”	**MS38-108** 6’
MS38-105 8”	**MS38-108** 8”
MS38-106 10”	**MS38-110** 10”
MS38-107 12”	**MS38-111** 12”



MS38-112 6” x8” **MS38-115** 6” x8”
MS38-113 8” x10” **MS38-116** 8” x10”
MS38-114 10” x12” **MS38-117** 10” x12”
With Lid Without Lid




MS38-118
Instruments Box
23x13x5cm




MS38-119
SCALER TRAY
10 x 20 mm



MS38-120
Dental PRF Box




Wire Mesh Cassettes
MS39-101-A (280x190x35mm) (20 instruments)
MS39-101-B (200x140x35mm) (10 instruments)
MS39-101-C (200x100x35mm) (07 instruments)
MS39-101-D (195x80x35mm) (05 instruments)




MS39-102-C
Dental Scalar Cassettes,A,B,C,D



MS39-103
Dental Scalar Clip Cassettes
for 10 Pcs



MS39-104
Dental Scalar Clip Cassettes
for 8 Pcs




MS39-105
Dental Scalar Cassettes,A,B,C,D




MS39-106
Dental Instruments Cassettes
for 5 Pcs



MS39-107
Dental Scalar Cassettes
for 10 Pcs



MS39-108
Dental Scalar Cassettes,A,B,C,D



The Oral Surgery Kit 1001:
1: 3.0mm Aspiration Tip / **2:** Serrated Lucas Surgical Curette, Small
3: Serrated Lucas Surgical Curette Large / **4:** Miller #1 Bone File / **5:** Cotton Pliers
6: Medium Straight Elevator / **7:** #77R Back Action Elevator **8:** Small Straight Elevato
9: Lower Universal Apical Forceps / **10:** Lower Cowhorn Forceps
11: Upper Universal Forceps / **12:** Lower Universal Forceps / 13: Cone Socket Mirror Handle
13: Bard-Parker Scalpel Blade Handle / **14:** Curved Kelly Haemostat / **15:** #5 Cone Socket Mirror
16: Crile-Wood Needle Holder / **17:** Seldin Flap Retractor / **18:** Molt 9 Periosteal Elevator
19: UNC Perio-Probe / #23 Explorer / **20:** Minnesota Retractor / **21:** Blumenthal Rongeur
22: Super-Cut Kelly Scissors / **23:** Adson/Brown Tissue Forceps, 2x8 Row Teeth
24: Extra Large Cassette



The Soft Tissue Kit # 2002
Description
1: Perio Chisel #13 / **2:** Large Cassette / **3:** Locking Cotton Pliers
4: Drop Blade Scalpel Handle / **5:** Cone Socket Mirror Handle
6: #5 Cone Socket Mirror / **7:** Molt Curette 2/4 / **8:** Crile-Wood Needle Holder / **9:** Sugarman 3/4 Perio File / **10:** Goldman-Fox #11 Perio Knife
11: Perio-Probe / #23 Explorer
12: Younger-Good 7/8 Curette
13: Goldman Fox Scissors



Micro Surgery Kit# 3003
1: Perio Chisel #13L / **2:** Perio Chisel/ **3:** Medium Cassette
4: Micro Cooley Forceps / **5:** Twist Joint Castro-Viejo Needle Holder
6: Periosteal Elevator / **7:** Castro-Viejo Scissors
8: Molt 9 Periosteal Elevator




The Implant Site Dev Kit # 4004
1: Small Distal Periotome TT Coated / **2:** Small Mesial Periotome TT Coated
3: Serrated Lucas Surgical Curette, Small / **4:** 3.5mm Graft Syringe
5: Large Cassette / **6:** Upper Universal Deep Grip Forceps
7: Lower Universal Deep Grip Forceps / **8:** Bone Condenser #12
9: Crile-Wood Needle Holder / **10:** Mini Buser Periosteal Elevator
11: Super-Cur Iris Scissors **12:** Suture Forceps, 1.6mm



Basic Implant Kit # 5005
Instruments Cassette
Crile Wood Needle Holders TC
La Grange Scissor TC
Adson Tissue Forceps TC
Periotomen Set
Scalpel Handle
Buser Periosteal Elevator
Prichard Periosteal Elevator
Lucas Bone Curette



MS Hygiene Kit 6006:
1: Medium Cassette / **2:** Cone Socket Mirror Handle
3: #5 Cone Socket Mirror / **4:** UNC Perio-Probe, Color
5: #204S Scaler / **6:** Gracey 11/12 Curette
7: Gracey 13/14 Curette / **8:** Columbia 4R/4L Curette
9: Gracey 1/2 Curette **10:** H6/7 Scale



Basic Exam Kit 7007
1: Small Cassette / **2:** Cotton Pliers
3: Cone Socket Mirror Handle
4: #5 Cone Socket Mirror
5: UNC Perio-Probe / #23 Explorer



Basic Exam Pack:8008
1: Regular #5 Mirror & Handle, Solid Handle
2: Explorer 17/23 , Solid Handle
3: Mariam Cotton Forceps

INSTRUMENTS CARE

Processing instruction for instruments

1STEP ONE:
Following clinical use, all instruments must be thoroughly rinsed with water to remove all residue prior to instrument processing.

STEP TWO:
Instruments should be placed into an ultrasonic bath for a period of 20 minutes to ensure complete instrument cleaning.
Note: Ultrasonic Solution should be changed daily or when the solution appears to be visibly soiled.

STEP THREE:
All instruments must be rinsed with water and completely dried before autociaving. All stainless-steel instruments should be placed in separate sterilization pouches and not to be mixed with any carbide instruments during instruments processing.

NOTE: Only demineralized water should be used in autoclave equipment according to autoclave manufacturer instructions

- ▶ All stainless-steel instruments must be cleaned and thoroughly dried prior to sterilization.
- ▶ We recommend using only a non-corrosive detergent.
- ▶ Abrasive brushes or materials should not be used to clean instruments.

2-STAINLESS STEEL INSTRUMENTS
We use many different types of stainless steels in the production of instruments dependent upon the design and application of the instrument. All stainless steels have a high nickel and chromium content to maximize corrosion resistance but will corrode and discolor when subjected to high concentrations of certain chemicals.

Do Not expose stainless steel dental instruments to the following chemicals. These chemicals will cause an adverse reaction and may destroy your instruments: Chlorine or Chlorinated products, Household Bleach, Tarter and Stain Remover, Aluminum Chloride, Aqua Regia, Barium Chloride, Bichloride of Mercury, Calcium Chloride, Carbolic Acid, Chlorinated Lime, Citric Acid, Dakin's Solution, Ferric Chloride, Ferrous Chloride, Hydrochloric Acid, Iodine, Lysol®, Mercury Chloride, Mercury Salts, Phenol, Potassium Permanganate, Potassium Thiocyanate, Sodium Hypochlorite (bleach), Stannous Chloride, Sulfuric Acid and Tartaric Acid (Tarter & Stain Remover).

Note: Be sure to only use solutions and chemicals that are compatible with stainless steel hand instruments.

Do Not use sodium chloride solution to clean the instruments. Never exceed temperatures 350°C F / 177°C C as this will have an adverse effect on the temper of the steel.

Recommended:
Steam sterilization in an autoclave is the benchmark process and is recommended. Sterilization requires a constant temperature of 134°C-134°CC for 5-10 minutes.

3 TITANIUM NITRIDE COATED INSTRUMENTS
Titanium coatings are used to increase the surface hardness of instrument tips and reduce "pullback" when manipulating composite materials. Proper care should be taken to remove any residual composite materials from the blade within 5 minutes after use by wiping the blade with alcohol on a 2×2 gauze. Composite materials or residue can harden on the blade and affect the quality of future restorations. Composite residue left to harden on the blade can not be removed without damaging the surface finish of the composite placement instrument. Titanium Nitride Coated instruments can be cleaned and sterilized using the same recognized acceptable methods as stainless steel instruments.

Note: Never use abrasives to clean the surface of any titanium coated composite placement instruments, elevators, periotomes or curettes.

HINGED INSTRUMENTS
All hinge instruments should be allowed to thoroughly dry in an open position after ultrasonic and manual cleaning procedures. We recommend that you only use lubricants such as surgical milk specifically formulated for dental and surgical instruments and follow manufacturer's instructions for applications.

Storage

Sterile, packaged instruments should be stored properly so that they are protected from dust, moisture, insects, vermin and extreme temperature.

Thank You.
MRSHAH CO. (GmbH)



Evolution Does not Imitate
It innovates